



Product Code	43116
Unit of sale	1 pair
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
--------	----

EAN/UPC	3701127751864
---------	---------------

PARAF	950079968
-------	-----------

Description

TIKI CRUTCHES - red/blue - pair

For children or teenagers. The handle adjusts from 16 to 21 cm.

Made in France.

Technical Specifications

Model TIKI

Weight kg: 0.35

Height handle to floor cm: 55-77

Material of handle: Polypropylene

Diameter upper/down pipe mm: 22/19

Number of adjustment of height: 10 times adjustable every 20 mm

Finishing/material: grey anodization/aluminium

Certified and tested for kg: 100

Standard used for testing: E N ISO 11334-1

Adjustable handgrip: 3 times

Handgrip size: Ø x l x h mm 32x100x160/210