

WALKER WITH 2 INTERNAL WHEELS



Product Code	43129
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2195550
CND	Y120606
EAN/UPC	8023279431292
UMDNS	14393
PARAF	977816420

Description

WALKING FRAME - 2 INTERNAL WHEELS

Made of lightweight aluminium tubing with padded handgrips, 2 non-slip ends and 2 with wheels. Foldable.
Adjustable height 80-90 cm.
Instructions: GB, FR, IT, SE.
Weight: 2.5 kg - Load: 100 kg.
Size 51x58 cm.