

MOVIGO WALKING FRAME - load 165 kg - anodised blue



Product Code	43133
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

Description

MOVIGO WALKING FRAME - load 165 kg - anodised blue - adult

This device provides support and stability for people who have difficulty walking.
MOVIGO is equipped with an inner reinforcement in the front foot area to strengthen the frame for maximum safety when walking.
Aluminium frame with polypropylene handles.
- adjustable height 81-91 cm
- lightweight
- reinforced tubes for enhanced rigidity
Max load: 165 kg.
Weight: 1.96 kg.

Supplied with multilingual manual: GB, FR, IT, ES, PT, DE, NL, PL, HU, SI, RO, TR.

Made in France.