



|               |                |
|---------------|----------------|
| Product Code  | 34314          |
| Unit of sale  | 1 pc           |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | I              |

|            |                |
|------------|----------------|
| UK-REP     | No             |
| CH-REP     | Yes            |
| RDM (NSIS) | 2837798        |
| CND        | Z12021099      |
| EAN/UPC    | 8023279343144  |
| PARAF      | 903147344      |
| GTIN       | 08023279343144 |
| GMDN       | 46827          |

## Description

Blade N° 3 MILLER - adults/large  
Overall Length: 195 mm  
Distal Width: 16.5 mm

- Great features of the unique GIMA blades and handles
- standard 2.5 V vacuum bulb
  - light ergonomic, satin finish handle to reduce glare
  - Stainless Steel contacts guarantee long working life
  - easy sterilization ETO or steam
  - full range of Mc Intosh & Miller blades, latex free
  - ISO 7376 fittings allow use of blades with existing handles
  - set is packaged in an anti-shock case and in cardboard box
  - 2.5 V handles work with both alkaline batteries or rechargeable batteries
  - good light transmission of 3,500 LUX with 2.5 V handles
  - autoclavable up to 134° for 5 minutes (approx. 2,000 times)
  - non magnetic blades