

# IDEAL ROLLATOR WITH 4 WHEELS



Product Code	43142
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes

RDM (NSIS)	2831155
CND	Y120606
EAN/UPC	8023279431421
UMDNS	17630
PARAF	977816471
GMDN	37951

## Description

IDEAL ROLLATOR - blue

Powder coated steel, adjustable height and foldable frame. With 4 solid castors (2 fixed rear and 2 swivel front). Ergonomic handgrip and hand brakes. Supplied with PU foam seat, reflector stickers, shopping basket and tray.

Weight: 10.8 kg - Load: 120 kg.  
Carton size: 590x255x715 mm.

Height adjustable 78-97 cm  
Tray 39 x 25 cm  
Wheels Ø 20 cm (8")  
Overall dimensions: 66 x 59.2 cm

Multilingual box and instructions: GB, FR, IT, ES