



Product Code	43153
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279431537
PARAF	950077154

Description

ELECTRIC PEDAL EXERCISER WITH DISPLAY

Made of ABS and steel with integrated LED display to track and read time, speed and calories. 6 adjustable speed levels.

Infrared remote control (battery not included).

Power supply: 80 W, 220V-50 Hz.

Multilingual box and manual: GB, FR, IT, ES, PT, DE, SE, LV, SI, GR, RO, PL.

Size: 39.5x34.5x26 cm

Packing size: 40x32x15 cm

Net weight: 4.2 kg

Non slip pedals with straps
 Forward/reverse to change pedal direction
 Easy to read monitor
 Non slid rubber pads

Compact exercise bikes for moderate cardiovascular workout. The regular use stimulates blood circulation and increases muscle strength. Can be placed on a tabletop to exercise arms or on the floor as a bicycle-type exerciser while sitting.