

FULLY ADJUSTABLE STANDUP ROLLATOR - aluminium - red



| | |
|---------------|----------------|
| Product Code | 43172 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|------------|---------------|
| RDM (NSIS) | 2329040 |
| CND | Y120606 |
| EAN/UPC | 8023279431728 |

Description

FULLY ADJUSTABLE STAND UP ROLLATOR

It helps adults to stand upright and look ahead with better support.

Compact one hand folding design, aluminum stable frame, dual hand brake system and fully adjustable armrest, backrest and length.

Detachable front wheels and stand up handles.

Padded seat and backrest for more comfort.

Two front (Ø 25.4 cm) and two rear (Ø 20.3 cm) PVC wheels, cane holder and shopping bag.

Packing size: 78x28.5xh 45 cm

Weight: 12 kg.

Load: 136 kg.

Dual hand brake system