



| | |
|---------------|----------------|
| Product Code | 43176 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|------------|--------------|
| UK-REP | No |
| CH-REP | Yes |
| RDM (NSIS) | 2955240 |
| CND | Y122406 |
| EAN/UPC | 805347072986 |

Description

UNIVERSAL BELTS FOR ARMCHAIRS, COMMODES, WHEELCHAIRS

- PREMIUM LINE ABDOMINAL BELT

Aids used for positioning and handling patients with safety.
Equipped with button acetal closure to increase safety level and prevent it from releasing.

Ideal for wheelchairs, commodes and armchairs (with 43178).
Coating in class 1 fire-retardant and water-repellent fabric.
Washable at max. 40°

Made in Italy.