



Product Code	43177
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes
RDM (NSIS)	2955225
CND	Y122406
EAN/UPC	8053470772955

## Description

UNIVERSAL BELTS FOR ARMCHAIRS, COMMODES, WHEELCHAIRS

- PREMIUM LINE HARNESS

Aids used for positioning and handling patients with safety.  
Equipped with button acetal closure to increase safety level and prevent it from releasing.

Ideal for wheelchairs, commodes and armchairs (with 43178).  
Coating in class 1 fire-retardant and water-repellent fabric.  
Washable at max. 40°

Made in Italy.