

TRANSAC TRANSFER BELT 120 cm - medium



Product Code	43182
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
--------	----

CND	Y123099
EAN/UPC	3701127727777
GMDN	041209

Description

TRANSAC TRANSFER BELT - medium

Polyurethane transfer belt to be placed around the patient's or the caregiver's waist, to help the patient sit or stand with minimum effort and maximum safety.

It has 5 vertical and 4 horizontal padded handles with polyester strap and nylon loop.

Length: 120 cm

Transferring load: 160 kg.

Washable at 60°.

Made in France.