

# RESTRAINT BELT L-XL 80-122 cm - single movement



Product Code	43185
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes

## Description

### BED RESTRAINT BELTS

- Single movement belt
- Size: L-XL
- Circumference cm: 80-122

100% cotton safety belts used for control of bedridden patient.  
Equipped with 1 safety lock that can be released with a special magnetic key.  
The belts allow a single movement to position the patient in a supine position.

Washable at 90°  
Belt height 21 cm.

Made in Italy.