



Product Code	43185
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2339614
CND	M030599
EAN/UPC	8053470771194

Description

BED RESTRAINT BELTS

- Single movement belt
- Size: L-XL
- Circumference cm: 80-122

100% cotton safety belts used for control of bedridden patient.
Equipped with 1 safety lock that can be released with a special magnetic key.
The belts allow a single movement to position the patient in a supine position.

Washable at 90°
Belt height 21 cm.

Made in Italy.