

# RESTRAINT BELT S-M 60-88 cm - double movement



Product Code	43186
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2339612
CND	M030599
EAN/UPC	8053470771156

## Description

### BED RESTRAINT BELTS

- Double movement belt
- Size: S-M
- Circumference cm: 60-88

100% cotton safety belts used for control of bedridden patient.

Equipped with 3 safety locks that can be released with a special magnetic key.

The belts allow a single or double movement to position the patient in a supine or lateral position.

Washable at 95°

Belt height 21 cm.

Made in Italy.