

RESTRAINT BELT L-XL 80-122 cm - double movement



Product Code	43187
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes

Description

BED RESTRAINT BELTS

- Double movement belt
- Size: L-XL
- Circumference cm: 80-122

100% cotton safety belts used for control of bedridden patient.
Equipped with 3 safety locks that can be released with a special magnetic key.
The belts allow a single or double movement to position the patient in a supine or lateral position.

Washable at 95°
Belt height 21 cm.

Made in Italy.