

COMFORT COMMODE CHAIR - height adjustable



Product Code	43191
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2044548
CND	Y091203
EAN/UPC	3401048584939
PARAF	977816521

Description

COMFORT COMMODE CHAIR - height adjustable - blue

Padded and removable seat for best comfort and high level of hygiene.
Pan is removable from both sides and from the top.

MAX LOAD 150 Kg

Made in France.

Technical Specifications

Total height: 81-96 cm
Armrests width: 53 cm
Total width: 61 cm
Height armrest-floor: 66-81 cm
Height seat-floor: 47-62 cm
Height seat-backrest: 34 cm
Seat plate depth: 44 cm
Total depth: 52 cm
Weight: 6.3 kg