



Product Code	43193
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
EAN/UPC	8023279431933

## Description

3 in 1 aluminium lifting chair which can be used as raised toilet seat, commode chair and shower chair. Designed with spring-loaded uplift mechanism and safety lock, that help users to stand up, pushing on the armrests. Self-powered, lean forward armrests for easy exit. Anti-slip tips with height adjustable legs. Delivered with bucket, bucket cover and splash-shield. Multilingual manual: GB, FR, IT, ES, PT, DE, SE, GR.

3 IN 1: COMMODE CHAIR, SHOWER CHAIR, TOILET RAISER

Cushion and backrest made of waterproof durable quality PU foam  
Adjustable seat height  
Self powered, lean forward armrest for easy exit  
Spring-loaded uplift mechanism  
Rugged, rust-proof aluminium frame

## Technical Specifications

Width: 52-55.8 cm  
Depth: 69.2-73 cm  
Height: 85-95.5 cm  
Packing: 64.5x58.5x39.5 cm  
Max load: 136 kg - Weight: 7.1 kg