

SMART COMMODE WHEELCHAIR



Product Code	43198
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	2370720
CND	Y122106
EAN/UPC	8023279431988

Description

Made of aluminium with anodized surface. Fixed height. Equipped with 4 castors with brakes, one large footrest, and pan with lid.
Indoor use. Multilingual manual: GB, FR, IT, ES, PT, DE, SE, PL, HR, RO, GR.

4 castors Ø 10 cm (4 ")

Packing size:
72x58x36 cm
Weight: 8.8 kg
Max load: 125 kg
Seat: 43x43 cm