

ABDOMINAL BELT for commode/wheelchair



Product Code	43216
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	1153308
CND	Y122406
EAN/UPC	8023279432169
PARAF	977816596

Description

UNIVERSAL BELTS FOR ARMCHAIRS, COMMODOES, WHEELCHAIRS

- SMART LINE ABDOMINAL BELT

Padded belt to prevent patient sitting on wheelchair from sliding forward.
Equipped with two black propylene webbing (height 5 cm), adjustable up to 100 cm.
One-hand release safety clip closure.
Harness has two adjustable suspenders with sturdy strap.
Washable at 30°.

Belt size: 48xh 15 cm, max. adjustable length: 248 cm.

Made in Italy.