



Product Code	43217
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	1153310
CND	Y122406
EAN/UPC	8023279432176
PARAF	977816608

Description

UNIVERSAL BELTS FOR ARMCHAIRS, COMMODES, WHEELCHAIRS

- SMART LINE HARNESS

Padded harness to prevent patient sitting on wheelchair from sliding forward.
Equipped with two black propylene webbing (height 5 cm), adjustable up to 100 cm.
One-hand release safety clip closure.
Harness has two adjustable suspenders with sturdy strap.
Washable at 30°.

Harness size: 60xh 37 cm, max. adjustable length: 260 cm.

Made in Italy.