



| | |
|---------------|----------------|
| Product Code | 43225 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|--------|----|
| UK-REP | No |
| CH-REP | No |

| | |
|------------|---------------|
| RDM (NSIS) | 2184716 |
| CND | Y122106 |
| EAN/UPC | 8023279432251 |
| UMDNS | 16949 |
| PARAF | 977816646 |

Description

TRANSIT WHEELCHAIR - 46 cm seat - black tissue

- double cross bar
- fixed armrest, PVC pad
- fixed footrest, flip-up plate, with leg strap
- 15 cm (6") solid front wheels
- 30 cm (12") pneumatic rear wheels
- locking brakes on rear wheels
- handles brakes for rear wheels
- folding back
- pocket on the back

Technical Specifications

- Total width: 63 cm
- Folded dimensions: -
- Seat width: 46 cm
- Rear wheel diameter: 53 cm
- Front wheel diameter: 15 cm
- Seat height: 49 cm
- Total wheelchair height: 91 cm
- Total wheelchair length: -
- Backrest height: -
- Seat depth: -
- Max load: 100 kg
- Kit size: 76 x 30 x 76 cm
- Net Weight: 14 Kg
- Gross Weight: 16.3 Kg