



Product Code	43235
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	2343875
CND	Y122103
EAN/UPC	8023279432350
UMDNS	16949

## Description

CLASSIC ALUMINIUM WHEELCHAIR - 46 cm seat - black tissue

- double cross bar
  - padded flipback armrest
  - detachable footrest, flip-up plastic plate with leg strap
  - 20 cm (8") solid PU front castors
  - 61 cm (24") solid PU rear wheels with quick release axles
- Multilingual manual: GB, FR, IT, ES, PT, SE.

Locking brakes on rear wheels.

## Technical Specifications

- Total width: 64 cm
- Folded dimensions: 37 cm
- Seat width: 46 cm
- Rear wheel diameter: 60 cm
- Front wheel diameter: 20 cm
- Seat height: 46 cm
- Total wheelchair height: 96 cm
- Seat depth: 42 cm
- Backrest height: 45.5 cm
- Max load: 120 kg
- Kit size: 82x31x96
- Net Weight: 13 kg
- Gross Weight: 14.5 kg