



Product Code	43250
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	2184715
CND	Y122106
EAN/UPC	8023279432503
UMDNS	16949
PARAF	977816774

## Description

QUEEN ALUMINIUM WHEELCHAIR - 46 cm seat - light blue tissue

- double cross bar
- extra light aluminium frame only 9.5 kg
- fixed armrest, cloth pad
- fixed flip-up footrest with leg strap
- 15 cm (6") front solid castors
- 30 cm (12") pneumatic rear wheels
- folding back
- locking brakes on rear wheels

Large pocket on the back

## Technical Specifications

- Total width: 60 cm
- Folded dimensions: 29.5 cm
- Seat width: 46 cm
- Rear wheel diameter: 30 cm
- Front wheel diameter: 15 cm
- Seat height: 48 cm
- Total wheelchair height: 90 cm
- Total wheelchair length: 100 cm
- Max load: 100 kg
- Kit size: 71 x 28 x 72 cm
- Net Weight: 9.5 Kg
- Gross Weight: 11 Kg