



| | |
|---------------|----------------|
| Product Code | 43261 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|------------|----------------|
| UK-REP | Yes |
| CH-REP | Yes |
| RDM (NSIS) | 2537327 |
| CND | Y122106 |
| EAN/UPC | 8023279432619 |
| UMDNS | 16949 |
| PARAF | 975866221 |
| GTIN | 08023279432619 |
| GMDN | 41512 |

Description

ESSEX WHEELCHAIR - 46 cm - black tissue

- blue powder coated steel frame
- rear solid PU wheels Ø 30.5 cm (12")
- front solid PU wheels Ø 20 cm (8")
- locking brakes on rear wheels
- front powder steel fork
- padded flip-up armrest
- detachable aluminium footrests with 6 height adjustable clips, flip-up plate, heel straps

Manuale: GB, FR, IT, ES, PT, DE, GR, SE

Large pocket on the back
Black nylon seat and back

Technical Specifications

- Total width: 61 cm
- Folded dimensions: 25 cm
- Seat width: 46 cm
- Rear wheel diameter: 30.5 cm
- Front wheel diameter: 20 cm
- Seat height: 50 cm
- Total wheelchair height: 91 cm
- Total wheelchair length: 96 cm
- Backrest height: 40 cm
- Seat depth: 41 cm
- Max load: 110 kg
- Kit size: 76 x 25 x 93 cm
- Net Weight: 16 Kg
- Gross Weight: 18 Kg