



Product Code	43262
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes
RDM (NSIS)	2537327
CND	Y122106
EAN/UPC	8023279432626
UMDNS	16949
PARAF	975866233
GTIN	08023279432626
GMDN	41512

Description

ESSEX WHEELCHAIR - 51 cm - black tissue

- blue powder coated steel frame
- rear solid PU wheels Ø 30.5 cm (12")
- front solid PU wheels Ø 20 cm (8")
- locking brakes on rear wheels
- front powder steel fork
- padded flip-up armrest
- detachable aluminium footrests with 6 height adjustable clips, flip-up plate, heel straps

Manual: GB, FR, IT, ES, PT, DE, GR, SE

Large pocket on the back
Black nylon seat and back
Double cross braces

Technical Specifications

- Total width: 66 cm
- Folded dimensions: 25 cm
- Seat width: 51 cm
- Rear wheel diameter: 30.5 cm
- Front wheel diameter: 20 cm
- Seat height: 50 cm
- Total wheelchair height: 91 cm
- Total wheelchair length: 96 cm
- Backrest height: 40 cm
- Seat depth: 41 cm
- Max load: 110 kg
- Kit size: 76 x 25 x 93 cm
- Net Weight: 16 Kg
- Gross Weight: 18 Kg