



|               |                |
|---------------|----------------|
| Product Code  | 43270          |
| Unit of sale  | 1 pc           |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | I              |

|        |     |
|--------|-----|
| UK-REP | Yes |
| CH-REP | Yes |

|            |               |
|------------|---------------|
| RDM (NSIS) | 2537327       |
| CND        | Y122106       |
| EAN/UPC    | 8023279432701 |
| UMDNS      | 16949         |
| PARAF      | 975866258     |
| GMDN       | 41512         |

## Description

OXFORD WHEELCHAIR - 43 cm - black tissue

- blue powder coated steel frame
- pneumatic rear wheels Ø 61 cm (24") with chromed steel handrim and quick release axles
- front solid PU wheels Ø 20 cm (8")
- locking brakes on rear wheels
- front powder steel fork
- padded flip-up armrests
- detachable footrests, flip-up plate, with leg straps

Black nylon padded seat and back  
Handgrip to fold wheelchair  
Quick release system

Supplied with pump and manual: GB, FR, IT, ES, PT, DE, GR, SE.

## Technical Specifications

- Total width: 60 cm
- Folded dimensions: 26.5 cm
- Seat width: 43 cm
- Rear wheel diameter: 61 cm
- Front wheel diameter: 20 cm
- Seat height: 50 cm
- Total wheelchair height: 91 cm
- Total wheelchair length: 103 cm
- Backrest height: 40 cm
- Seat depth: 44 cm
- Max load: 120 kg
- Kit size: 82 x 27 x 93 cm
- Net Weight: 18 kg
- Gross Weight: 20 Kg