



Product Code	43274
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2537327
CND	Y122106
EAN/UPC	8023279432749
UMDNS	16949
PARAF	975866284
GTIN	08023279432749
GMDN	41512

## Description

OXFORD PLUS WHEELCHAIR - 46 cm - black tissue - with solid PU rear wheels Ø 61 cm (24") and anti tipper system

- blue powder coated steel frame
- solid PU rear wheels Ø 61 cm (24") with chromed steel handrim and quick release axles
- front solid PU wheels Ø 20 cm (8")
- locking brakes on rear wheels
- front powder steel fork
- padded flip-up armrests
- detachable footrests, flip-up plate, with leg straps

Black nylon padded seat and back  
Handgrip to fold wheelchair  
Quick release system  
Large pocket on the back  
Anti tipper system

Supplied with pump and manual: GB, FR, IT, ES, PT, DE, GR, SE.

## Technical Specifications

- Total width: 63 cm
- Folded dimensions: 26.5 cm
- Seat width: 46 cm
- Rear wheel diameter: 61 cm
- Front wheel diameter: 20 cm
- Seat height: 50 cm
- Total wheelchair height: 91 cm
- Total wheelchair length: 103 cm
- Backrest height: 40 cm
- Seat depth: 44 cm
- Max load: 120 kg
- Kit size: 82 x 27 x 93 cm
- Net Weight: 18 kg
- Gross Weight: 20 Kg