



Product Code	43275
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2537327
CND	Y122106
EAN/UPC	8023279432756
UMDNS	16949
PARAF	975866296
GTIN	08023279432756
GMDN	41512

Description

OXFORD PLUS WHEELCHAIR - 51 cm - black tissue - with solid PU rear wheels Ø 61 cm (24") and anti tipper system

- blue powder coated steel frame
- solid PU rear wheels Ø 61 cm (24") with chromed steel handrim and quick release axles
- front solid PU wheels Ø 20 cm (8")
- locking brakes on rear wheels
- front powder steel fork
- padded flip-up armrests
- detachable footrests, flip-up plate, with leg straps

Double cross braces
Black nylon padded seat and back
Handgrip to fold wheelchair
Quick release system
Large pocket on the back
Anti tipper system

Supplied with pump and manual: GB, FR, IT, ES, PT, DE, GR, SE.

Technical Specifications

- Total width: 68 cm
- Folded dimensions: 26.5 cm
- Seat width: 51 cm
- Rear wheel diameter: 61 cm
- Front wheel diameter: 20 cm
- Seat height: 50 cm
- Total wheelchair height: 91 cm
- Total wheelchair length: 103 cm
- Backrest height: 40 cm
- Seat depth: 44 cm
- Max load: 120 kg
- Kit size: 82 x 27 x 93 cm
- Net Weight: 18 kg
- Gross Weight: 20 Kg