



Product Code	43283
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	2324847
CND	Y122103
EAN/UPC	8023279432831

Description

FOLDABLE COMMODE WHEELCHAIR

Chrome steel structure device ideal for patients with low mobility. It has fixed armrests and footrests, solid castors and sturdy rear wheels. U shape commode seat is extremely durable, suitable for indoor use. Supplied with toilet bowl.

Rear wheels diameter: 58 cm (22").
Front castors diameter: 20 cm (7").

Technical Specifications

- Total width: 64 cm
- Seat width: 44 cm
- Rear wheel diameter: 58 cm
- Front wheel diameter: 20 cm
- Seat height: 50 cm
- Total wheelchair height: 87 cm
- Total wheelchair length: 105 cm
- Backrest height: 37 cm
- Max load: 100 kg
- Kit size: 97 x 29 x 89 cm
- Net Weight: 19 kg
- Gross Weight: 21,3 Kg