

# MAXLITE F.O. PEDIATRIC SET 3 BLADES MILLER 00-0-1 - 2.5 V - 8.000 lux



Product Code	34332
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	283778
CND	Z12021003
EAN/UPC	8023279343328
PARAF	903147433
GTIN	08023279343328
GMDN	15076

## Description

### MAXLITE F.O. XENON/HALOGEN SET - 2.5 V

- Pediatric set (3 blades Miller 00-0-1)

- F.O. laryngoscopes provide cool bright illumination
- MaxliTE F.O. blades with integrated fibre
- Xenon/Halogen F.O. rechargeable battery handle
- 2.5 V rechargeable battery handle (battery not included). 2.5 V rechargeable battery handle can be used with standard or rechargeable batteries. Rechargeable battery handles should be used with charging station 31542 (not included)
- non magnetic blades