

MAXLITE F.O. ADULT SET 3 BLADES MC-INTOSH 1-2-3 - 2.5 V - 8.000 lux



Product Code	34333
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	283776
CND	Z12021003
EAN/UPC	8023279343335
PARAF	903147458
GTIN	08023279343335
GMDN	15076

Description

MAXLITE F.O. XENON/HALOGEN SET - 2.5 V

- Adult set (3 blades Mc Intosh 1-2-3)

- F.O. laryngoscopes provide cool bright illumination
- MaxliTE F.O. blades with integrated fibre
- Xenon/Halogen F.O. rechargeable battery handle
- 2.5 V rechargeable battery handle (battery not included). 2.5 V rechargeable battery handle can be used with standard or rechargeable batteries. Rechargeable battery handles should be used with charging station 31542 (not included)
- non magnetic blades