

MAXLITE F.O.LED PEDIATRIC SET 3 BLADES MILLER 00-0-1 - 2.5 V - 20.000 lux



Product Code	34336
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	283778
CND	Z12021003
EAN/UPC	8023279343366
GTIN	08023279343366
GMDN	15076

Description

MAXLITE F.O. LED SET - 2.5 V

- Pediatric set (3 blades Miller 00-0-1)

- F.O. laryngoscopes provide cool bright illumination
- Maxlite F.O. blades with integrated fibre
- LED F.O. rechargeable battery handle
- 2.5 V rechargeable battery handle (battery not included). 2.5 V rechargeable battery handle can be used with standard or rechargeable batteries. Rechargeable battery handles should be used with charging station 31542 (not included)
- non magnetic blades