



| | |
|---------------|----------------|
| Product Code | 34340 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | II A |

| | |
|--------|-----|
| UK-REP | Yes |
| CH-REP | Yes |

| | |
|------------|----------------|
| RDM (NSIS) | 1831045 |
| CND | Z1203020408 |
| EAN/UPC | 8023279343403 |
| PARAF | 903147460 |
| GTIN | 08023279343403 |
| GMDN | 45607 |

Description

Small and lightweight wrist Oximeter available for home and clinical long-term and comfortable monitoring. It can be simply worn on the patient's wrist as a watch to start the monitoring and can be used for sports to measure SpO2 and pulse rate (PR). Not suitable for pediatric use.

- audible and visual alarm with adjustable alarm limits
- trend review function of SpO2 and PR
- selectable storage modes: 1 group/1s, 1 group/4s, and 1 group/8s
- up to 500 hours storage for SpO2 and PR in default mode (1 group/8s)
- clock setting and date display function
- integrated software to download and store all measurements on PC
- removable finger probe
- comfortable to be worn on wrist as a watch