

HYDRAULIC PATIENT TRANSFER CHAIR



Product Code	43430
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	2692857
CND	V080502
EAN/UPC	8023279434309
PARAF	950850800

Description

HYDRAULIC PATIENT TRANSFER CHAIR

Hydraulic transfer chair made of steel, ideal for transferring the patient from bed to toilet, shower or other sites.

It includes 3 interchangeable seat panels:

- fixed toilet panel that can be used with or without a removable toilet pan
- removable plastic shower panel
- removable soft cushion seat panel

Shower and soft cushion seat panels can be disassembled whenever necessary to accommodate varied conditions, leaving the area more hygienic and cleaner.

- 180° easy opening and closing, flexible for movement
- hydraulic lifting adjustment up to 15 cm
- load capacity: 125 kg
- 4 directional castors (2 frontal with brake)

Packing size: 79x63x33 cm.
Seat size: 48x37xh 42-57cm

Multilingual manual: GB, FR, IT, ES, PT, DE, SE, PL, HU, HR, RO, GR.