

HYDRAULIC PATIENT LIFTER - load 150 kg



Product Code	43450
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2517667
CND	V0805030101
EAN/UPC	8023279434507

Description

HYDRAULIC PATIENT LIFTER - max load 150 kg

Hydraulic lifter made of steel frame, surface with coating.
Front Ø 7.5 cm wheels without brake, rear Ø 10 cm wheels with brake.
The legs can be opened by pedal.
Height adjustment by hydraulic pressure or manual control.
Supplied with sling (43464).

Multilingual manual: GB, FR, IT, ES, PT, SE, HR, RO.

Technical Specifications

Overall width: 590 mm
Hook height: 950-1650 mm
Push handle height: 1140 mm
Base open width: 860 mm
Base clearance: 115 mm
Base length: 1190 mm
Arm length: 990 mm
Load: 150 kg