

ELECTRIC PATIENT LIFTER - load 150 kg



Product Code	43455
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
--------	-----

RDM (NSIS)	2517666
CND	V0805030102
EAN/UPC	8023279434552

Description

ELECTRIC FOLDABLE PATIENT LIFTER - max load 150 kg

Electric lifter made of steel frame, surface with coating.
Foldable for easy storage and transfer.
Front Ø 7.6 cm wheels with brake, rear Ø 10 cm wheels with brake.
Legs can be opened by pedal.
Hand controller with one function.
It has emergency stop button, emergency drop device and backup battery.

Supplied with sling (43464).

Multilingual manual: GB, FR, IT, ES, PT, SE, HR, RO.

Technical Specifications

Overall length: 1200 mm
Overall width: 580 mm
Foldable size: 1250x580x310 mm
Boom height: 0-1730 mm
Push handle height: 1130 mm
Base open width: 830 mm
Base clearance: 115 mm
Base length: 1120 mm
Base inner width: 475-730 mm
Power supply: 220V - 50Hz
Load: 150 kg