



Product Code	43460
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	2517665
CND	V0805030102
EAN/UPC	8023279434606

Description

ELECTRIC ALUMINIUM PATIENT LIFTER - FOLDABLE

Electric lifter with aluminum frame and coating surface.
 Foldable for easy storage and transfer.
 7.5 cm Ø wheels: front wheels without brake, rear wheels with brake.
 The legs can be opened by pedal.
 One function hand controller.
 Provided with emergency stop button, emergency drop device and backup battery.
 Supplied with sling (43464).

Multilingual manual: GB, FR, IT, ES, PT, SE, HR, RO.

Technical Specifications

Overall length: 1180 mm
 Overall width: 710 mm
 Foldable size: 1310x710xh 335 mm
 Boom height: 0-1610 mm
 Push handle height: 1160 mm
 Base open width: 1070 mm
 Base clearance: 100 mm
 Base length: 1120 mm
 Base inner width: 600-960 mm
 Power supply: 220 V - 50 Hz
 Load: 150 kg