



Product Code	43500
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes

RDM (NSIS)	2537393
CND	V0299
EAN/UPC	8023279435009
GMDN	38140

Description

NEONATAL CRADLE with trolley

Transparent unbreakable plastic cradle with round edges and ventilation holes.
Epoxy powder painted mainframe made of curved and welded tubular steel.
Trolley shape allows horizontal and Trendelenburg inclinations and bedside positioning.
4 synthetic rubber castors Ø 75 mm, two with brake.
Weight: 15 kg.
Max load: 10 kg.
Manual: GB, FR, IT, ES.
Delivered disassembled.

Dimensions: 83.5x52.5xh 90 cm

Made in Italy.