



Product Code	43502
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	No
CH-REP	No
EAN/UPC	8023279435023

## Description

### BABY CRADLE - light blue

Cradles with stylish and ergonomic design to provide a comfortable experience for babies.

Transparent baby basin made of plastic.

Wear-resistant, antiskid, shock absorption double side castors Ø 76 mm (3") can freely rotate 360°.

Easy to operate handle to control gas spring.

12° reverse Trendelenburg.

Adjustable height by hydraulic pump.

Size: 875x525xh 650 to 900 mm.

Packing size: 930x560x850 mm.

Supplied with baby basin.