

# CHILDREN BED - 2 cranks - light blue



Product Code	43506
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2321126
CND	V080699
EAN/UPC	8023279435061

## Description

### CHILDREN BED - 2 cranks - light blue

Beds for children with frame made of sheet metal bending and welding.  
Plastic parts made of ABS material and epoxy painting.  
Bed ends are made of PP blow moulding to protect bed and patient.

#### Guardrail

Aluminium alloy liftable guardrail with 8 cm distance among protective tubes. 3 gears to adjust height, anti-pinch hands to protect children.

#### IV Pole

Stainless steel IV pole, bearing 15 kg

#### Hand control

Manual cranks to control back and leg liftings

#### Platform

Concave vent hole design of the panel allows strong air permeability and can effectively prevent bedsores

#### Shoes' holder

Utility and shoes holder

#### Castors Ø 12.7 cm (5")

Castors with brakes

## Technical Specifications

Size: 1880x900xh 1260 mm  
Bed height: 585 mm  
Mattress platform: 1770x720 mm  
Back-rest adjustment: 0-75° (± 10°)  
Knee-rest adjustment: 0-35° (± 10°)

## Standard accessories

1 set full length side rail  
4 buffer air springs  
1 stainless steel IV pole bearing 15 kg  
Shoe holder  
4 IV pole provision (2 on each head)

