

CHILDREN BED - 2 cranks - pink



Product Code	43507
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2321126
CND	V080699
EAN/UPC	8023279435078

Description

CHILDREN BED - 2 cranks - pink

Beds for children with frame made of sheet metal bending and welding.
Plastic parts made of ABS material and epoxy painting.
Bed ends are made of PP blow moulding to protect bed and patient.

Guardrail

Aluminium alloy liftable guardrail with 8 cm distance among protective tubes. 3 gears to adjust height, anti-pinch hands to protect children.

IV Pole

Stainless steel IV pole, bearing 15 kg

Hand control

Manual cranks to control back and leg liftings

Platform

Concave vent hole design of the panel allows strong air permeability and can effectively prevent bedsores

Shoes' holder

Utility and shoes holder

Castors Ø 12.7 cm (5")

Castors with brakes

Technical Specifications

Size: 1880x900xh 1260 mm
Bed height: 585 mm
Mattress platform: 1770x720 mm
Back-rest adjustment: 0-75° (± 10°)
Knee-rest adjustment: 0-35° (± 10°)

Standard accessories

1 set full length side rail
4 buffer air springs
1 stainless steel IV pole bearing 15 kg
Shoe holder
4 IV pole provision (2 on each head)

