

DOUBLE BASE UNIT 4 drawers + 1 door with 2 shelves - white



Product Code	43510
--------------	-------

Unit of sale	1 pc
--------------	------

Minimum order	1
---------------	---

Type	No medical device
------	-------------------

Class	
-------	--

UK-REP	No
--------	----

CH-REP	No
--------	----

EAN/UPC	8023279435108
---------	---------------

Description

Double base units

- 4 drawers + door with 2 shelves - white

High quality double base units made of wood covered with scratch-resistant and anti-reflective melamine finish.
The top, bottom and back of the structure and the shelves are covered with abs edges.
4 legs (height: 15 cm) with levellers.

Two models: one with 4 drawers and door with 2 shelves (41.4 cm) and the other one with 2 doors and 2 shelves (86.4 cm).
It can be requested in different colours.

Made in Italy.