

# DOUBLE BASE UNIT 2 doors with 2 shelves - white



Product Code	43515
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	No
CH-REP	No
EAN/UPC	8023279435153

## Description

### Double base units

- 2 doors with 2 shelves - white

High quality double base units made of wood covered with scratch-resistant and anti-reflective melamine finish. The top, bottom and back of the structure and the shelves are covered with abs edges. 4 legs (height: 15 cm) with levellers.

Two models: one with 4 drawers and door with 2 shelves (41.4 cm) and the other one with 2 doors and 2 shelves (86.4 cm). It can be requested in different colours.

Made in Italy.