

# ROOM WARDROBE - streaked beige - disassembled



Product Code	43525
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	No
CH-REP	No

EAN/UPC	8023279435252
---------	---------------

## Description

### ROOM WARDROBE - streaked beige

Structure made of melamine.

Internal central partition that forms two compartments: one compartment with 4 adjustable shelves and drawer, other compartment has an upper shelf with rod for clothes and bottom shelf to store shoes.

Provided with cylinder lock, 4 visible feet (height: 15 cm) with levellers.

Size: 90x54xh 185 cm.

Delivered disassembled.

Made in Italy.