

OLYMPUS BEDSIDE TABLE with TILTING OVERBED TABLE



Product Code	43552
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

EAN/UPC 8023279435528

Description

OLYMPUS BEDSIDE TABLE with TILTING OVERBED TABLE

Double-sided bedside tables made of polyethylene.
The inner compartment has a removable plastic shelf and a removable grid at the bottom to guarantee ventilation.
Side panels include plastic towel holder, two cup and two bottle holders all equipped with retaining clamps.
Drawer and doors made of moulded polyethylene plastic material with integrated handles.
Height adjustable overbed tables with epoxy powder-coated steel frame with raised edges to avoid spills.
360° rotation and horizontal maneuverability of plate.
Tilting top at different angles.
Four 55 mm Ø twin swivel castors.

Made in Italy.

Technical Specifications

Size: 61.5x48.5xh 91 cm
Weight: 21.5 kg
Drawer compartments size: small: 36x10xh 8.5 cm / large: 36x18xh 8.5 cm
Interior compartment size: 32.5x38xh 44 cm

Overbed table:
size: 57x44 cm
height: from 91 to 121 cm
horizontal extension: 15 cm