

# HEINE STANDARD F.O. HANDLE 3.5V with recharg. Battery



Product Code	34372
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
RDM (NSIS)	1379511
CND	R9099
EAN/UPC	8023279343724
UMDNS	15-076

## Description

### HEINE RECHARGEABLE LARYNGOSCOPE HANDLES AND BATTERIES

- Standard F.O. handle 3.5 V Li-ion with rechargeable battery
- HEINE Code: F-007.22.894

Need NT4 charger unit (code 34374).

The "Classic+" range of F.O. blades with HEINE HiLite®

Fibre bundle with over 6,500 individual micro-fibres offer improved light transmission and longer life. Together with the innovative HEINE LEDHQ light in the handle or the traditional XHL Xenon Halogen Bulb, the exclusive HEINE HiLite® Fibre Optic bundle with up to 4.3 mm of Ø offers 3X or 40% brighter light as with other conventional halogen bulbs.

5-years guarantee or 4,000 standard autoclave cycles of 134°C/5 min. with more than 1,000 lux on all blades with HEINE HiLite® fibre bundle.

#### Long lasting bulbs

HEINE XHL® Xenon Halogen bulbs incorporate the latest Xenon Halogen Technology to give you up to 40% more light! XHL 3.5 V bulbs last up to twice as long as 2.5 V bulbs and give up to three times the light of conventional 2.5 V Halogen laryngoscope bulbs. HEINE LEDHQ handles offer even three times the light intensity of the 3.5V Xenon Halogen bulbs for an outstanding and constant bright illumination.

#### Long lasting 3.5 V rechargeable batteries

HEINE 3.5 V systems (LEDHQ or XHL) use a long-lasting 3.5 V lithium rechargeable battery. HEINE F.O. laryngoscope blades and handles are compatible with the ISO 7376 (green standard).