

# SPLINTER WITH MAGNIFYING GLASS 2X - 8 cm



Product Code	24391
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes

RDM (NSIS)	2835530
CND	L99
EAN/UPC	8023279243918
PARAF	977810643
GTIN	08023279243918

## Description

SPLINTER WITH MAGNIFYING GLASS 2X - 8 cm

Magnifying glass splinter tweezer made of stainless steel.