

2-SECTION WOODEN MASSAGE TABLE - black



Product Code	44000
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279440003
PARAF	977816812

Description

2-SECTION WOODEN MASSAGE TABLE - black

Two-section foldable height adjustable massage table with German beech wood frame and polyurethane cover with 5 cm foam thickness. Tables are supplied with headrest, hanging armrest, mouth-nose plug and carrying bag.

Product open size: 185x70xh 62/87 cm
Packing size: 95x20xh 75 cm
Net/Gross weight: 14/17 kg
Maximum working load: 227 kg (500 lbs)
Manual: GB, FR, IT, ES, PT, DE