

# 2-SECTION WOODEN MASSAGE TABLE - cream



Product Code	44003
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required

EAN/UPC	8023279440034
PARAF	977816848

## Description

### 2-SECTION WOODEN MASSAGE TABLE - cream

Two-section foldable height adjustable massage table with German beech wood frame and polyurethane cover with 5 cm foam thickness. Tables are supplied with headrest, hanged armrest, mouth-nose plug and carrying bag.

Product open size: 185x70xh 62/87 cm  
Packing size: 95x20xh 75 cm  
Net/Gross weight: 14/17 kg  
Maximum working load: 227 kg (500 lbs)  
Manual: GB, FR, IT, ES, PT, DE