

HEINE STANDARD HANDLE 3.5V LED Li-ion



Product Code	34402
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
--------	-----

EAN/UPC	4053755129928
UMDNS	15-076

Description

F-008.22.891 Standard LED handle 3.5 V Li-ion with recharg. battery (Need NT300 charger unit code 34449)

We introduce the „+“ range of F.O. blades with HEINE HiLite® Fiber Optic bundle with over 6.000 fibre optic and the state of the art HEINE XHL® Xenon Halogen bulbs in the Handle. The innovative HEINE HiLite® Fiber Optic bundle with up to 100 % more cross-section improves light transmission and makes the blade more resistant to autoclaving than ever. HEINE offers a 5-year guarantee or 4000 standard autoclave cycles of 134°C/5 min. with more than 1000 lux on all blades with HEINE HiLite® fibers. HEINE XHL® Xenon Halogen bulbs incorporate the latest Xenon Halogen Technology to give you up to 40 % more light! All HEINE F.O.+ laryngoscope systems offer exceptional value for money and effectively save money for you. HEINE XHL® Xenon Halogen 3.5 V bulbs last up to twice as long as 2.5 V bulbs and give up to three times the light of conventional 2.5 V Halogen laryngoscopes bulbs. HEINE 3.5 V systems use a long-lasting 3.5 V rechargeable battery which avoids the need to replace non-rechargeable batteries at frequent intervals. HEINE F.O. laryngoscope blades and handles are compatible with the ISO 7376 (green standard).