

# 2-SECTION ALUMINIUM MASSAGE TABLE - black



Product Code	44020
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279440201
PARAF	977816875

## Description

### 2-SECTION ALUMINIUM MASSAGE TABLE - black

Two-section foldable, height adjustable massage table with aluminium frame and polyurethane cover with 5 cm thickness. Tables are supplied with headrest, side armrests, hanged armrest, mouth-nose plug and carrying bag.

Product open size: 186x70xh 59/80 cm  
Packing size: 95x20xh 75 cm  
Net/Gross weight: 13/16 kg  
Maximum working load: 180 kg (397 lbs)  
Manual: GB, FR, IT, ES, PT, DE